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info@bionexia.de · www.bionexia.de

Rev. 01/01-2007

Quick test for qualitative detection of human hemoglobin in stools

General information · application and evaluation

Application

FOB plus is a visual immunochromatographic quick test for the qualitative detection of human hemoglobin in stool samples. The test is only intended for professional in vitro diagnostic use. The test was designed to supply evidence of illnesses in the large intestine, such as bowel cancer and larger hemorrhoids.

Bowel cancer is one of the most frequently diagnosed types of cancer and one of the most frequent causes of death from cancer (Lieberman, 1994; MMWP, 1995).

Bowel cancer can be detected early by means of tests for hidden blood in the stools and the mortality rate can therefore be reduced. (Dam et. al., 1995; Miller, 1995; and Lang, 1996).

Previous tests for hidden blood in stools used the Guajak method that requires a special diet to avoid false positive and negative results. The highly specific FOB plus has therefore been developed to detect human hemoglobin in stool samples. The test is based on an immunochemical method that improves on the detection of bowel illnesses in comparison with the traditional Guajak test by means of higher sensitivity and specificity (Frommer et. al., 1988; St. John et. al., 1993).

Test principle

The immunological Bionexia® FOB quick test is based on a immunochromatographic method, in which only human hemoglobin is recognized through specific antibody reactions.

The hemoglobin contained in the stool sample reacts with specific monoclonal antibodies that are attached to gold particles. This complex is distributed on the membrane and reaches the test band (T) on which the anti-hemoglobin is present. In the case of a positive result, the molecules from the stool sample loaded with gold-marked antibodies attach to the test band (T) and are visible by means of a pink/red coloration. In the case of a negative result, there are no hemoglobin molecules that can attach to the test band (T) as complexes and therefore, there can also be no coloration of the test band (T).

If the control strip (C) turns red/pink in color, this is a guarantee that the sample has been taken correctly and has migrated correctly. The test is therefore valid.

Storage and stability

The test can be stored at 2° - 30° in a sealed bag until the expiration date.

Please note

- Valid only for in vitro diagnostic use and only for professional use.
- Do not use the test after the expiration date.
Do not mix test collection containers with different batches.
- Only open the test packaging if the test is actually going to be performed.
- All patient samples should be handled as if they could carry disease.
- The extraction solution contains a small amount of sodium azide.
- Patients should follow the instructions for obtaining a sample exactly and should not take any sample during menstruation or in case of bleeding hemorrhoids.
Use of the enclosed stool container is recommended in order to avoid mixing with blood from urine samples or interfering components in the toilet water.

Materials supplied:

- 25 individually packed test cassettes
- 25 sample-taking tubes
- 25 sets of instructions on taking stool samples
- 25 stool containers
- 1 set of usage instruction

Sample taking and pre-treatment


1. Store the tube at room temperature before use.
2. Collect the stool using the stool container.
3. Hold the tube vertically and open the sealing cap. Remove the sample-taker.
4. Place the spiral rod at three different locations in the stool sample.
5. Place the sample-taker with the stool sample into the tube and seal well.
6. Place the tightly-sealed tube in the small bag.
7. Fill out the cover sheet from these short instructions with the name and date of birth of the patient.
8. Store the tube in as cool a place as possible overnight and protected from light.
9. Stool samples in the buffer solution are stable at room temperature for up to 15 days. Within this period of time, it is possible to store, transport and test the sample at room temperature without needing to take further measures. For later tests, the sample must be stored at +2° to +8°C for the entire period of time and be brought to room temperature shortly before the test is performed.
10. If the sample is taken at the patient's home, the sample should be stored there protected from light and in a cool place.

Performance of the test

Quality control

1. The test contains a procedural control. A colored line that appears in the control region shows that each test was performed correctly.


It is easy for the background to become slightly yellowish in color during testing, depending on the color of the stool sample. This is acceptable, as long as the evaluations of the test results are not adversely affected.

Procedure


1. The test cassette and the patient's sample should be brought to room temperature (20°C to 30°C) before testing.
2. First take the test cassette out of the packaging when you perform the test. Bring the cassette to room temperature to avoid moisture deposits in the membrane. Mark the cassette with the name of the patient or with another form of identification.
3. Shake the stool sample tube carefully, in order to ensure that the stool sample mixes properly with the extraction solution.
4. Unscrew the protective cap. Take a paper towel and break the seal of the collection container with the rotary motion. Hold the collection container upright and add 4 drops of the solution into the round sample window of the cassette.
5. Read the results after 5 minutes. Strongly positive results may be evaluated even sooner. Evaluate the results within 8 minutes after applying the Solution on the round sample window.

Evaluation of the results

Positive: 2 pink-colored lines appear; one in the control region (C) and one in the test region (T).

If you are testing strongly positive samples, the intensity of the control line may be reduced. It is not recommended to compare the intensity of the lines.

Negative: Only 1 pink-colored line appears in the control region (C). No line appears in the test region (T).

Invalid: If no red line appears in either region, this is a sign that the test is not functioning properly, or that the test materials are not correct. In this case, repeat the test with a new test cassette or contact the manufacturer for technical support.

Scope of application

1. The test is used for the qualitative determination of human hemoglobin in stool samples. A positive result occurs if human hemoglobin is present in the sample. In addition to intestinal bleeding, blood in stools may also have other causes, such as hemorrhoids or the mixing of the stool sample with urine containing blood.
2. Not all intestinal bleeding are the result of benign or malignant polyps. The information that is obtained by this test should be examined in conjunction with other clinical test methods.
3. Negative results cannot exclude the possibility of bowel cancer, because some polyps and types of bowel cancer are intermittent or do not bleed at all. In addition, blood is not evenly distributed in stool samples. Bowel polyps in the early stages do not necessarily bleed.
4. Urine and excessive dilution of samples with water from the toilet bowl can lead to false test results. It is therefore recommended that the stool container be used.
5. Stool samples should not be taken during menstruation or 3 days before or after, in the case of bleeding caused by constipation, bleeding hemorrhoids, or in the case of medicine being administered rectally. This could lead to false positive results.

This test is less sensitive in the detection of bleeding in the small intestine because a resorption of the blood can take place during further process through the bowel.

Analytical capacity

A. Analytical sensitivity

A sample, that contains human hemoglobin at a concentration of approx. 2 µg/g stool (= 40 ng hemoglobin/mL buffer after extraction), generates a positive result. In some cases, samples that contain human hemoglobin at a concentration of less than 4 ng /mL in extraction buffer), can still test positive. This is particularly true after duration of 30 minutes. Prozone effect Bionexia FOBplus also works reliably at extremely high hemoglobin values (>500,000ng/mL).

B. Specificity

Bionexia® FOB plus is specifically for human hemoglobin and has no cross-reactivity with hemoglobin from cattle, pigs, rabbits, horses, or sheep at a concentration of up to 0.5 mg/mL in the extraction buffer.

FOB plus also displays no cross-reactivity with Bilirubin, vitamin C, and horseradish peroxidases.

Note

A study recently published by the Shinshu University School of Medicine in Japan investigated the cost/benefit ratio of multiple tests. The study showed that the relative sensitivity increases with an increased number of tests, although the relative specificity decreases.

Results:

Number of tests	Sensitivity	Specificity
1	58%	96%
2	89%	95%
3	100%	94%

Literature:

1. Dam, J.V., et. al.; Fecal Blood Screening for Colorectal Cancer; Archive of Internal Medicine; (1995) 155:2389-2402
2. Frommer, D.J. et. al.; Improved Screening for Colorectal Cancer by Immunological Detection of Occult Blood; British Medical Journal; (1988) 296:1092-1094
3. Lieberman, D.; Screening/Early Detection Model for Colorectal Cancer, Why Screen? Cancer Supplement; (1994) 74(7) : 2023-2027
4. Miller, A.B.; An Epidemiological Perspective on Cancer Screening; Clinical Biochemistry (1995) 28(1) : 41-48
5. Ransohoff, D.F. and Lang, C.A.; Improving the Fecal Occult- Blood Test: The New England Journal of Medicine; (1996) 334 (3) : 189-190
6. Screening for Colorectal Cancer-United States, 1992-1993, and New Guidelines; Mobility and Mortality Weekly Report; (1995) 45 (5) : 107-110
7. St. John, D.J.B., et. al.; Evaluation of New Occult Blood Test for Detection of Colorectal Neoplasia; Gastroenterology; (1993) 104:1661-1668
8. Yamamoto M.; Nakama H.; Cost-effectiveness analysis of immunochemical occult blood screening for colorectal cancer among three fecal sampling methods : Hepatogastroenterology; 2000 Mar-Apr; 47(32):396-9.

Explanation of symbols

LOT = batch

 = Expiration date

Do you have questions on the use or principle of the test?
Contact your supplier or the manufacturer.

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